

# Configure Your Wheelchair

## Step 1: Select Chair Type

### Pedal

Wheelchairs that enable an individual to safely move and be independent.

### Positioning

Tilt-in-space plus recline functionality for the highest level of comfort and care.

### Hygiene

Standard

Bariatric

Revive Shower Commode

Revive Bariatric Shower Commode

### Glider

Adult

Pediatric

Tranquille Glider

Aspire Glider

### Rehab

Durable and versatile wheelchairs for users in rehab and therapy.

Flip to other side

## Step 2a: Select Positioning Wheelchair (Caregiver-Propelled)

Features, Functionality, and Customization Level

Bariatric

Vanguard Bariatric Wheelchair



Basic

Best



Centric Positioning Wheelchair



Midline Positioning Wheelchair



Elite Positioning Wheelchair Elite Tilt Wheelchair



Synthesis Positioning Wheelchair Synthesis Transport Wheelchair

250 lbs. weight capacity  
Up to 25° Tilt  
Up to 45° Recline

350 lbs. weight capacity  
Up to 38° Tilt  
Up to 90° Recline

350 lbs. weight capacity  
Up to 40° Tilt

350 lbs. weight capacity  
Up to 40° Tilt  
Up to 90° Recline

### Optional Add-ons:

Additional Positioning Padding (APP)

Huntington's Special Padding (HSP)

WC-19 Transport Package

Acute Care Configuration (ACC)

Slim Posture Padding (SPP)

## Step 2b: Select Pedal Wheelchair (Self-Propelled [Foot or Hand])

Features, Functionality, and Customization Level

Basic

Best



Sashay Pedal Wheelchair



Latitude Pedal Wheelchair with Dynamic Rocking Feature



Encore Pedal Wheelchair



Encore Pedal Wheelchair with Dynamic Rocking feature

350 lbs. weight capacity  
Up to 18° Tilt

250 lbs. weight capacity  
Up to 33° Tilt

350 lbs. weight capacity  
Up to 25° Tilt

250 lbs. weight capacity  
Up to 33° Tilt

### Optional Add-ons:

Slim Posture Padding (SPP)

Additional Positioning Padding (APP)

## Step 3: Select Accessories, Colors, & Customizations

Visit our website or call us for accessories, colors, and customization options.  
[www.brodaseating.com](http://www.brodaseating.com) 1-844-552-7632

# Configure Your Rehab Chair

More Independent

Less Independent

## Step 2c: Select Chair Frame

**BRODA®**



*Encore Rehab Wheelchair*



*Latitude Rehab Wheelchair*



*Comfort Rehab Wheelchair*



*Elite Rehab Wheelchair*

250 lbs. weight capacity Up to 33° Tilt Dynamic Rocking HCPCS E1161	250 lbs. weight capacity Up to 33° Tilt Dynamic Rocking HCPCS E1161	250 lbs. weight capacity Up to 45° Tilt HCPCS E1161	350 lbs. weight capacity Up to 40° Tilt Up to 80° Recline HCPCS E1161
--	--	---	--

For convenience, all Matrx® options ship on chair.

## Step 2c.ii: Select Seating

More Pressure Reduction

Less Pressure Reduction

**matrx®**  
SEATING SERIES



*Libra*



*Vi*



*PSP*



*PS*

HCPCS E2624/E2625	HCPCS E2607/E2608	HCPCS E2607	HCPCS E2605/E2606
-------------------	-------------------	-------------	-------------------

**Encore Rehab or Comfort Rehab**

**Latitude Rehab or Elite Rehab**

## Step 2c.iii: Select Back

More Supportive

More Accommodating

**matrx®**  
SEATING SERIES



*Elite TR*




*Elite & Elite Deep*



*PB & PB Deep*

HCPCS E2620/E2621	HCPCS E2620/E2621	HCPCS E2615	HCPCS E2601/E2611
-------------------	-------------------	-------------	-------------------

**BRODA®**



*Comfort Tension Seating®  
Additional Positioning Padding (APP)  
Slim Posture Padding (SPP)*

## Step 2c.iv: Select Headrest

**matrx®**  
SEATING SERIES



*Elan Headrest*

HCPCS E0955/E0128
-------------------

## Step 3: Select Additional Accessories

<b>Seating &amp; Back Options</b>	Pelvic Obliquity Kit, Cushion Rigidizer, Seat Wedges, Abductors, Fluid Supplements, Fixed and Swing Away Laterals, Privacy Flap
<b>Chair Accessories</b>	Swing Away, Removable, and Elevating Legrests, Removable and Swing Away Armrests

*Call us for more options!  
1-844-552-7632*

HCPCS codes supplied are for reference only and do not ensure coverage or payment. The HCPCS provided are based on our interpretation of the code definitions. HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.